



QL[®]

INFLATION

DEVICES

Atrion[®]
Medical

CUSTOM BUILT

The QL[®] Inflation Device

Perfect for a wide array of applications such as balloon catheterization, stent deployment and fluid delivery.

Available in a variety of configurations so you can be sure you get exactly what you need.



Smart Design

Stable design and ergonomic handle make it easy to control. The quick latch mechanism allows for rapid deflation, even from max pressure. Clear barrel provides medical staff with a precise view of the content and volume allowing user to perform thorough purging.

Broad Applications

- Cardiology
- Peripheral
- Radiology
- Urological
- Gastrointestinal
- Neurology
- Orthopedics
- Biliary
- Pulmonary
- ENT

Atrion[®]
Medical

FOR YOUR APPLICATION

QL[®] Inflation Device at a Glance

The QL[®] is available in two types



Quick Latch



Press-to-Release

- 510(k) Cleared
- CE Marked
- Ergonomic Handle
- Clear Barrel
- Available Sterile or Non-Sterile
- Made in the U.S.A.

Stock Configurations

Part Number	Size	Pressure	Type
96500	10cc	15 ATM	Quick Latch
96501	10cc	30 ATM	Quick Latch
96351	14cc	55 ATM	Quick Latch
96502	20cc	30 ATM	Quick Latch
96450	20cc	30 ATM	Press-to-Release
96358	25cc	30 ATM	Quick Latch
96311	25cc	40 ATM	Quick Latch
96313	40cc	15 ATM	Quick Latch
96314	60cc	15 ATM	Quick Latch

Custom Ordering

Available in a variety of configurations. Custom name and logo can be printed on the gauge, carton and label. Minimum order quantities apply.

Unique

If your application has special requirements, Atrion will work with your team to develop a unique, custom inflation device, locking syringe or fluid delivery device to meet your needs.

Service

We work with our customers to ensure orders are built to specified requirements and delivered on time.

Regulatory Approvals

Approved for use in U.S., Canada, and Europe.

+1 256-317-2123

info@atrionmedical.com

1426 Curt Francis Road, Arab, AL 35016

STANDARD
CUSTOM
UNIQUE



Atrion Medical Products, Inc. | 1426 Curt Francis Road, Arab, AL 35016

+1 (256) 317-2123 | info@atrionmedical.com | www.atrionmedical.com

Protected by one or more of the following U.S. Patents: 6,796,959 • 6,938,319 • 8,191,457 • 8,499,681.
Including Foreign Patents: CA 2,437,056, CN Z1 2009-10211532.9, CN Z1 2012-10429157.7, DE 60 2002 226 968.7,
DE 60 2009 010 578.7, DE 60 2008 030 736.0, EP 1385569, EP 2027846, EP 2586485, JP 4040469, JP 4233587,
JP 5117953, JP 5909768, JP 5909773, KR 10-152636, and Patents Pending.
QL® is a registered trademark of Atrion Medical Products, Inc.