



A Nordson MEDICAL Company



QL®

Inflation Devices & Fluid Dispensing Syringes

Atrion, now a part of Nordson MEDICAL, offers a versatile portfolio of high-pressure Inflation Devices and Fluid Dispensing Syringes. Our devices are engineered for a broad range of clinical applications including balloon catheter dilation, stent deployment, and high pressure fluid delivery. These devices are available in multiple configurations to support diverse procedural needs, ensuring clinicians have the tools they need to optimize outcomes. All products are manufactured in an ISO 13485, FDA registered facility and can be customized or private labeled to meet specific user requirements.

FEATURES & BENEFITS

● Smart Design

Designed for stability and ease of use, the ergonomic handle ensures precise control during operation. The quick-latch mechanism enables rapid deflation, even from maximum pressure levels. A transparent barrel offers clear visibility of contents and volume, allowing medical personnel to perform thorough purging with confidence.

● Custom-Engineered Solutions

Atrion's engineering team collaborates closely with customers to develop specialized inflation devices, locking syringes, and fluid delivery systems tailored to unique clinical or procedural requirements. Whether your application demands custom functionality, materials, or performance characteristics, Atrion offers flexible design and manufacturing capabilities to deliver solutions that align precisely with your specifications.

Custom Branding Options

- This product is available in multiple configurations to accommodate a variety of clinical requirements. Custom branding options are also offered, allowing you to print your name and logo on the gauge, carton, and label.

● Customer-Focused Service You Can Count On

Atrion partners closely with customers to ensure every order is built to exact specifications and delivered on time. Our commitment to precision and reliability helps support your success from production to delivery.

● Broad Applications

- Cardiology
- Peripheral Vascular
- Structural Heart
- Urology
- Gastroenterology
- Neurology
- Orthopedics
- Cell Delivery
- Pulmonary
- ENT

QL® Inflation Devices

Stock Configurations

Part Number	Size	Pressure	Type
96500	10cc	15 ATM	Quick Latch
96501	10cc	30 ATM	Quick Latch
96351	14cc	55 ATM	Quick Latch
96502	20cc	30 ATM	Quick Latch
96450	20cc	30 ATM	Press-to-Release
96358	25cc	30 ATM	Quick Latch
96311	25cc	40 ATM	Quick Latch
96313	40cc	15 ATM	Quick Latch
96518	45cc	15 ATM	Quick Latch
96314	60cc	15 ATM	Quick Latch
96522	60cc	20 ATM	Push Button

Devices for Balloon Catheter Dilations

QL1455 (14cc/55ATM)

Ultra High-Pressure



QL2540 (25cc/40ATM)

Ultra High-Pressure



QL45(45cc/15ATM)

Reduced Deployment Force



QL6015(60cc/15ATM)

Large Capacity



QL6020PB(60cc/20ATM)

Large Capacity

